


Apollo^{NC}

Stronger • Better • Safer

Non-Compliant, High Pressure
Coronary Balloon Catheter



High Pressure •
Non-Compliant •
6F Guide Kissing •
3.0% Axial Growth •
0.55% Compliance •

Apollo is a non-compliant high pressure post-dilatation catheter offering numerous features for enhanced and robust performance.

- Exceptionally minimal axial growth during inflation
- 5F guiding catheter compatible
- Versatile; Optimal balloon diameters and lengths to meet clinical demands (2.0-5.0 mm diameters) (8-18 mm lengths)
- Average balloon burst pressure is 30 atm

Brasmed
embrace life

Balloon Specification

Balloon Diameter (mm)	2.0, 2.25, 2.5, 2.75, 3.0, 3.25, 3.5, 3.75, 4.0, 4.5, 5.0
Balloon Length (mm)	8, 10, 12, 15, 18
Nominal Pressure	12atm
Rated Burst Pressure	22atm (ϕ 2.0-4.0mm) 20atm (ϕ 4.5-5.0mm)
Average Burst Pressure	30atm
Compliance	Non-compliant, 0.55%
Balloon Pleats	3
Tip Length (mm)	2.0/2.5
Distal Shaft OD	2.55F/2.6F
Proximal Shaft OD	2.0F
Guiding Catheter Compatibility	5F Compatible/6F compatible for kissing
Working Length (cm)	140
Coating	Slider™ Hydro from tip to adjacent guide wire port (external surface); Eel™ coating for guide wire lumen

Ordering Information

Balloon Diameter	Balloon Length/Product Code				
	8mm	10mm	12mm	15mm	18mm
2.0mm	802-2008	802-2010	802-2012	802-2015	802-2018
2.25mm	802-2208	802-2210	802-2212	802-2215	802-2218
2.5mm	802-2508	802-2510	802-2512	802-2515	802-2518
2.75mm	802-2708	802-2710	802-2712	802-2715	802-2718
3.0mm	802-3008	802-3010	802-3012	802-3015	802-3018
3.25mm	802-3208	802-3210	802-3212	802-3215	802-3218
3.5mm	802-3508	802-3510	802-3512	802-3515	802-3518
3.75mm	802-3708	802-3710	802-3712	802-3715	802-3718
4.0mm	802-4008	802-4010	802-4012	802-4015	802-4018
4.5mm	802-4508	802-4510	802-4512	802-4515	802-4518
5.0mm	802-5008	802-5010	802-5012	802-5015	802-5018



BrosMed Medical Co., Ltd.

15th Building, SMES Venture Park, Songshan Lake Hi-Tech Industrial Development Zone, Dongguan523808, Guangdong, China
Tel: +86 (769) 2289 2018 Fax: +86 (769) 2289 2016
E-mail: sales@brosmed.com

BrosMed Medical B.V.

Mgr. Buckxstraat 8, 6134 AP Sittard, The Netherlands
Office: +31(0) 850 140 901
Mail: cs@brosmed.com

FDA 510K



GRA-6002(CE)/Rev.09