


Conger Hydrophilic Guide Wires



PRODUCT FEATURES

- One-piece nitinol core construction promotes pushability and easy advancement
- One to one torque for navigating distal and angulated branch vessels
- Durable hydrophilic coating to facilitate crossing of tight lesions
- Stiff and standard versions accommodate a wide range of clinical applications
- Configurations: 45, 80, 150, 180, 230, 260 cm – 0.025", 0.032", 0.035", 0.038"

ORDERING INFORMATION

Catalogue No.	Tip Shape	Hardness	OD(inch)	Length(cm)
HNS-025-45	Straight	Standard	0.025	45
HNA-025-45	Angled	Standard	0.025	45
HNS-025-80	Straight	Standard	0.025	80
HNA-025-80	Angled	Standard	0.025	80
HNS-025-150	Straight	Standard	0.025	150
HNA-025-150	Angled	Standard	0.025	150
HNS-025-180	Straight	Standard	0.025	180
HNA-025-180	Angled	Standard	0.025	180
HNS-025-260	Straight	Standard	0.025	260
HNA-025-260	Angled	Standard	0.025	260
HNS-032-150	Straight	Standard	0.032	150
HNA-032-150	Angled	Standard	0.032	150
HNS-032-180	Straight	Standard	0.032	180

Catalogue No.	Tip Shape	Hardness	OD(inch)	Length(cm)
HNA-032-180	Angled	Standard	0.032	180
HNS-032-230	Straight	Standard	0.032	230
HNA-032-230	Angled	Standard	0.032	230
HNS-032-260	Straight	Standard	0.032	260
HNA-032-260	Angled	Standard	0.032	260
HNS-035-45	Straight	Standard	0.035	45
HNA-035-45	Angled	Standard	0.035	45
HNS-035-80	Straight	Standard	0.035	80
HNA-035-80	Angled	Standard	0.035	80
HNS-035-150	Straight	Standard	0.035	150
HNA-035-150	Angled	Standard	0.035	150
HXS-035-150	Straight	Stiff	0.035	150
HXA-035-150	Angled	Stiff	0.035	150
HNS-035-180	Straight	Standard	0.035	180
HNA-035-180	Angled	Standard	0.035	180
HXS-035-180	Straight	Stiff	0.035	180
HXA-035-180	Angled	Stiff	0.035	180
HNS-035-230	Straight	Standard	0.035	230
HNA-035-230	Angled	Standard	0.035	230
HXS-035-230	Straight	Stiff	0.035	230
HXA-035-230	Angled	Stiff	0.035	230
HNS-035-260	Straight	Standard	0.035	260
HNA-035-260	Angled	Standard	0.035	260
HXS-035-260	Straight	Stiff	0.035	260
HXA-035-260	Angled	Stiff	0.035	260
HNS-038-150	Straight	Standard	0.038	150
HNA-038-150	Angled	Standard	0.038	150
HNS-038-180	Straight	Standard	0.038	180
HNA-038-180	Angled	Standard	0.038	180