

# Connecting Tubings



## INTRODUCTION

The product is intended to connect pipelines and devices which provides pathway of the liquid-injection or/and invasive blood pressure monitoring in interventional operation.

## PRODUCT FEATURES

- High pressure resistance: 500psi/1200psi
- Transparent material makes the flow visible
- Customized sizes are available

## NON-BRAIDED TUBING

**500 PSI** female to male luer lock

External OD (mm)	Length (cm)						
	10	20	30	60	90	120	150
1.8	BM-CT1.80-10	BM-CT1.80-20	BM-CT1.80-30	BM-CT1.80-60	BM-CT1.80-90	BM-CT1.80-120	BM-CT1.80-150
2.0	BM-CT2.00-10	BM-CT2.00-20	BM-CT2.00-30	BM-CT2.00-60	BM-CT2.00-90	BM-CT2.00-120	BM-CT2.00-150
2.7	BM-CT2.70-10	BM-CT2.70-20	BM-CT2.70-30	BM-CT2.70-60	BM-CT2.70-90	BM-CT2.70-120	BM-CT2.70-150
3.0	BM-CT3.00-10	BM-CT3.00-20	BM-CT3.00-30	BM-CT3.00-60	BM-CT3.00-90	BM-CT3.00-120	BM-CT3.00-150
3.2	BM-CT3.20-10	BM-CT3.20-20	BM-CT3.20-30	BM-CT3.20-60	BM-CT3.20-90	BM-CT3.20-120	BM-CT3.20-150
3.65	BM-CT3.65-10	BM-CT3.65-20	BM-CT3.65-30	BM-CT3.65-60	BM-CT3.65-90	BM-CT3.65-120	BM-CT3.65-150
4.0	BM-CT4.00-10	BM-CT4.00-20	BM-CT4.00-30	BM-CT4.00-60	BM-CT4.00-90	BM-CT4.00-120	BM-CT4.00-150

## BRAIDED TUBING

**1200 PSI** female luer lock to male rotator

External OD (mm)	Length (cm)				
	20	60	90	120	150
3.65	BM-CT3.65-20B	BM-CT3.65-60B	BM-CT3.65-90B	BM-CT3.65-120B	BM-CT3.65-150B

## NON-BRAIDED TUBING

**1200 PSI** female luer lock to male rotator

External OD (mm)	Length (cm)				
	20	60	90	120	150
3.20	BM-CT3.20-20A	BM-CT3.20-60A	BM-CT3.20-90A	BM-CT3.20-120A	BM-CT3.20-150A