

# Freepath Hydrophilic Sheath Introducer Sets



## PRODUCT FEATURES

- Smooth transition between dilator and sheath to avoid puncture trauma
- Superior hydrophilic coating on both dilator and sheath
- High kink resistance
- Sophisticated tapered tip for easy entry
- Available lengths: 7, 10, 12, 25 cm
- Color coded : 4F—Red      5F—Grey  
                                6F—Green      7F—Orange

## DESIGNED FOR RADIAL ACCESS

Catalogue No.	Size(Fr.)	Length(cm)	Nitinol guide wire	I.V. Cannula	Syringe
H-N-4T	4F	12cm	0.025" × 45cm	20G × 5cm	2.5ml
H-N-5T	5F				
H-N-6T	6F				
H-N-7T	7F				

\*CE not available

## DESIGNED FOR FEMORAL ACCESS

Catalogue No.	Size(Fr.)	Length(cm)	Nitinol guide wire	I.V. Cannula	Syringe				
H-S-47	4F	7cm	0.025" × 45cm	18G × 6.5cm	2.5ml				
H-S-57	5F								
H-S-67	6F								
H-S-77	7F								
H-S-510	5F	10cm	0.035" × 45cm			18G × 6.5cm	2.5ml		
H-S-610	6F								
H-S-710	7F								
H-N-4	4F	12cm	0.025" × 45cm					18G × 6.5cm	2.5ml
H-N-5	5F								
H-N-6	6F								
H-N-7	7F								
H-L-5	5F	25cm	0.035" × 70cm						
H-L-6	6F								
H-L-7	7F								

\*CE not available