

Manifolds



PRODUCT FEATURES

- Compatible with most infusion pumps, syringes and intravenous infusion tubes to meet complex clinical needs
- Excellent pressure resistance and sealing performance
- Leak free and smooth handling
- On / Off types are available for different user habits

ONE PORT

Catalogue No.	Product description/Pressure
BM-MF1-ON	“ON” 500psi
BM-MF1-OFF	“OFF” 500psi

TWO PORTS

Catalogue No.	Product description/Pressure
BM-MF2-ON	“ON” 500psi
BM-MF2-OFF	“OFF” 500psi

THREE PORTS

Catalogue No.	Product description/Pressure
BM-MF3-ON	“ON” 500psi
BM-MF3-OFF	“OFF” 500psi