

Radial Artery Tourniquets



BM-RAT2



BM-RAT1

INTRODUCTION

The device is used to achieve hemostasis after the radial arterial catheterization.

PRODUCT FEATURES

- Unique design for radial artery hemostasis post percutaneous puncture
- Transparent fixed plates ensures clear visual control
- Soft and comfortable adjustable band

PRODUCT COMPOSITION

The device consists of a fixed plate, adjustable band and elastomer

| Catalogue No. | Product description |
|---------------|---------------------|
| BM-RAT1 | Cord button |
| BM-RAT2 | Rotary |